





#### **Tortillas – A Universal Pleasure**




- Increasing demand for convenient food
- May be used as a bread replacement
- Low in carbohydrates
- Rich in minerals, vitamins, and proteins
- Many possibilities, adapted to specific diet (GF, vegan,...) or requirements (organic, clean label)
- People are looking to increase the diversity of their meal
- Example: Pizza is the favorite dish in the world, Tortilla arrives at Position 10 after Croissant and Ramen Source: Taste Atlas

Bosc-Bierne, TIA Europe Conference 2023

## **Baking is a Transformation Process**




## **Potential Influence of Raw Material Quality (Flour)**



#### WATER ABSORPTION

capacity must be controlled to reach the right CONSISTENCY without excessive STICKINESS.

The dough must be sufficiently EXTENSIBLE so that it does not break, and NOT VERY ELASTIC so that it does not shrink after cutting.



The final product must have the good SHAPE, THICKNESS and COLOR. Conform ROLLABILITY and FOLDABILITY are required.

No DEFECT should be visible: No CRACKS nor PILLOWING

+ SHELF LIFE is a key parameter.

Bosc-Bierne, TIA Europe Conference 2023



> There are many different types of tortillas and great variation in manufacturing processes!

Wheat flour is the major and most-relevant ingredient used for manufacturing flour tortillas and the quality of the finished product depends greatly on the quality of the flour.

Wheat flour requirements are determined by the desired tortilla characteristics, the formula, processing conditions, and equipment.

Different types of flours are required for different tortilla processes; hot-press, die-cut, and handstretch. Dough preparation and ingredient utilization vary among different operations. Each operation involves a unique dough-forming procedure that then requires specific flour characteristics.

Bejosano and Alviola, 2015





Selecting the right flour requires careful consideration because *protein* and *starch properties* that cause longer shelf stability are the same factors that make tortilla diameter smaller.

A strong-protein flour makes tortillas with smaller diameter but with longer shelf stability. On the other hand, a weak protein-strength flour makes larger-diameter tortillas that have short shelf stability. Thus, flour with intermediate protein quality would be appropriate.

Both protein content and strength of flour show negative correlations with tortilla diameter.

Nevertheless, the quality of the end product also depends on other processing variables and the formulation that the manufacturer uses.

Bejosano and Alviola, 2015






Starch damage in flour affects the properties of tortillas.

Wheat flours with a low damaged starch content are reported to produce tortillas with a larger diameter and a better texture than those with a high damaged starch content.

There is an **ideal level of** damaged starch according to the level of proteins and the ideal area for accepting incoming flours should be defined.





#### **Starch damage**

 When starch damage occurs at low levels (or at optimum/moderate levels for the desired product) of the total wheat starch, they are considered to be beneficial to flour performance;

because they are readily accessible to water and enzyme penetration and serve the dual function of increasing flour absorption and providing fermentable maltose for the yeast.



#### Starch damage

#### **EFFECTS ON WATER ABSORPTION**

- Protein absorbs 1.8 times its weight of water,
- Pentosans: 10 times,
- Native starch: 0.4 times,


Damaged starch: 3 to 4 times!

#### MILLERS' OBJECTIVE

- Produce the largest amount of flour from a lot of wheat (economic profitability)
- Produce flours according to the specifications of their customers (high hydration capacity flours for example for bread making flours)

#### **BAKERS' INTEREST**

Add max content of water (better yield and better profit)







#### **Starch damage**

#### **EFFECTS ON DOUGH PROCESSING**

- Doughs produced with high mechanical starch modification levels may not be fully developed at Peak time (optimum mixing time) and display a second peak – Modify mixing time?
- Possibility of having a sticky dough probably because of weakened protein filler interfacial interactions and larger voids within the protein network?

Hackenberg, et al., 2019

 Damaged starch increases initial water absorption and prevents optimum gluten formation during mixing!!!



## Why do we measure starch damage?

#### **Starch damage**

#### **EFFECTS ON STALING/SHELF LIFE OF THE BREAD**

The **content of damaged starch** also influences the bread staling through increasing the retrogradation of the amylopectin and crumb firmness, and therefore affects <u>negatively</u> the quality of the resulting fresh products and their <u>shelf life</u>.

Le'on, et al., 2006





Effect of damaged starch on wheat flour tortilla quality




#### What is damaged starch?


# Structurally disruption/modification of starch granules in wheat and flour!


Starch damage, damaged starch, mechanically activated starch, mechanical starch modification

Small, Round (B-starch)



(Wilson, et al., 2006) Large, Lenticular (A-starch)







#### Where does it come from?

There are two types of starch damage:

 Damaged starch comes from the wheat itself (as a function of hardness) - genetic criteria –

The first is that which results in the starch granule being broken in two. Although the granule is clearly damaged, this type of damage results in starch that is not susceptible to attack by fungal  $\alpha$ -amylase.






#### Where does it come from?

Damaged starch also comes from the milling process - mechanical criteria -

The second more classic starch damage produced during milling results in granules that have partially or completely lost their crystallinity and are susceptible to fungal  $\alpha$ -amylase.



Barrera, et al., 2013





#### Starch damage during milling

 Roller milling damages a small but significant number of the starch granules in the flour.

In general, the level of damage depends not only on the hardness of the wheat used but also on the settings of the milling rolls. Indeed, given a constant wheat supply, flour with increasing water-absorption characteristics can be produced by adjusting the roll settings such that increasing levels of starch damage are produced.




Soft



### Starch damage during milling

## Impact of roll adjustments

Starch damage = f (Roll gap)






### Starch damage during milling

### Impact of roll adjustments

Damaged starch = f (Rolls pressure)








## **Objective**

- To assess the impact of mechanically damaged starch on processing and quality attributes of flour tortillas!
  - Hard and soft wheat





#### **Materials and Methods**





## **Materials and Methods: Damaged Starch Determination**


# AACC 76-33.01 Damaged Starch — Amperometric Method by SDmatic

 Based on the measurement of iodine absorption, it works on 1 gram of flour and provides results in only 10 minutes.



- Preparation of an iodine solution.
- Iodine creates an electrical flow.
- When the iodine fixes on the damaged starch, the intensity of the current decreases.
- 4. The less intense the electrical current, the higher the damaged starch content.








# **Materials and Methods: Tortilla Formula and Processing**


Ingredients	% f.b.
Flour	100 (14%MC)
Distilled water	Abs.
Shortening	7
Sugar	0.5
Salt	1.5
Sodium bicarbonate	1
SAP	1.8
SSL	0.5
Potassium sorbate	0.4
Calcium propionate	0.5

Liu, Hou, Cardin, Marquart and Dubat, 2017





# **Tortilla Processing - Hot-Press Method**





# **Materials and Methods: Tortilla Quality Evaluation**



Diameter

Thickness

Weight

Color (Lab)

Rollability

Texture






# **Results: Damaged Starch, Water Absorption**

	Soft Wheat Flour		Hard Wheat Flour		
	As is	1 Pass	As is	1+2 Pass	6+7 Pass
Damaged Starch Content, % (AACC 76-33.01)	4.62 ±0.15 d	5.83 ±0.16 c	12.46 ±0.06 b	12.87 ±0.24 b	13.70 ±0.08 a
Water Absorption, %	53.5	54.9	62.6	63.7	62.1




### **Results: Tortilla Dimensions**






# **Results: Rollability Day 1 versus Day 7**



✓ Effect of damaged starch on wheat flour tortilla quality



## Results: Resistance/Toughness Day 1 versus Day 7





Hard wheat flour tortilla: Increased starch damage led to increased toughness at Day 1. This difference was completely lost after 7 days of storage.

Soft wheat flour tortilla: Damaged starch content caused no difference in toughness at both test days, but at 7<sup>th</sup> day.

Toughness is the maximum resistance to probe penetration and extensibility is measured from distance prior to burst.



## **Results: Extensibility Day 1 versus Day 7**

At Day 1 there was no difference in extensibility between all the flours tested.

At Day 7, the tortillas produced with soft wheat flour and the hard wheat flour at 1+2 passes had the highest had extensibility.



Toughness is the maximum resistance to probe penetration and extensibility is measured from distance prior to burst.







#### **Discussion & Conclusion**

- > Tortilla diameter is negatively correlated with flour strength and protein content.
- Tortillas made from flours with more damaged starch had a smaller diameter but longer shelf stability (Arora 2003).

Wheat flours with a low damaged starch content are reported to produce tortillas with a lager diameter and a better texture than those with high damaged starch content (Wang et al., 2020).



#### **Discussion & Conclusion**

- ➤ Both toughness and extensibility are related to starch retrogradation and starchprotein interactions.
- The lower damaged starch content and polyphenol oxidase activity in Eastern US soft winter wheat, ESW when compared with hard wheat flour would be advantageous because ESW wheat would require less water for dough preparation and would produce larger and brighter-colored tortillas (Ma and Baik, 2023).





### **Materials**

#### 15 wheat varieties have been selected


		SKCS	Flour	
#	Variety	Hardness Index	Moisture Content	Flour Protein content
1	ALONDRA F2014	59,1	13,17	10,6
2	BACOREHUIS F2015	61,5	13,33	12,1
3	BORLAUG100 F2014	64,0	13,29	11,9
4	CIANO M2018	68,7	12,83	12,0
5	CISNE F2014	64,2	12,75	10,4
6	CONATRIGO F2015	63,0	12,79	11,7
7	FUERTEMAYO F2016	70,1	12,87	12,0
8	HANS F2019	66,5	12,9	11,7
9	KRONSTAD F2014	71,5	13,11	12,6
10	LUMINARIA	65,5	12,9	12,3
11	NORESTE F2018	66,2	13,14	12,6
12	NORMAN	58,3	12,71	11,8
13	TACUPETO F2001	62,6	12,84	11,5
14	VALLES F2015	76,3	12,63	12,5
15	VILLA JUAREZ	64,2	12,9	10,9
	Minimum	58,3	12,6	10,4
	Maximum	76,3	13,3	12,6
	Average	65,4	12,9	11,8



- All hard wheat varieties, with strong or medium strong gluten
- Grown under optimal conditions at the CIMMYT research station in the North of Mexico



#### **Methods**











Alveolab Standard test at constant hydration (ISO 27971)





Mixolab Chopin+ test (ISO 17718)





Tortilla Production
Test protocol
(internal method
– see next slide)






## **Tortilla Production Test Protocol (details)**

 A SIMPLE FORMULATION AND A STRAIGHT FORWARD TEST PROTOCOL





A COMPLETE TORTILLA ANALYSIS







Effect of damaged starch on wheat flour tortilla quality

## **Overview of the Tortilla Results**










Tortilla Quality Parameters								
#	WEIGHT (g)	DIAMETER (cm)	THICKNESS (mm)	OPACITY (%)	ROLLABILITY	FOLDABILITY	TOTAL SCORE	CONCLUSION
1	45,4	21,7	1,4	95	5	1	2	BAD
2	43,7	23,9	1,5	85	5	4	7	GOOD
3	43,6	22,9	1,5	85	5	3	4	MEDIUM
4	42,1	25,2	1,6	50	5	3	5	MEDIUM
5	43,5	23,5	1,6	80	4	2	5	MEDIUM
6	43	23,6	1,5	80	4	4	6	MEDIUM
7	43,2	22,3	1,4	90	4	2	3	BAD
8	44	22,8	1,4	75	4	2	1	BAD
9	43,7	23,9	1,3	85	4	3	6	MEDIUM
10	43,9	22,4	1,5	85	3	3	2	BAD
11	45,8	22,6	1,7	80	4	2	4	MEDIUM
12	45,1	23,2	1,5	95	5	4	7	GOOD
13	45	22,9	1,6	95	4	5	6	GOOD
14	45,7	21,7	1,5	90	4	3	6	MEDIUM
15								
Min	42,1	21,7	1,3	50	3	1	1	
Max	45,8	25,2	1,7	95	5	5	7	




## Damaged Starch, the Best Indicator for Overall Quality of Wheat for Tortilla Making



- •UCD values > 22 predictive of poor quality tortillas!
- •UCD values < 19 predictive of high quality tortillas!









The SDmatic allows the definition of simple specifications to select high quality wheat for making tortillas!

## **Study Conclusion**

### **OVERALL QUALITY**





**TORTILLA** 

**QUALITY** 

"Starch damage is the most important trait influencing the overall tortilla quality and should always be analyzed when producing flour for tortilla making" Vega and Ibba, TIA Europe 2021

**ROLLABILITY + OVERALL QUALITY** 











- While the tortilla quality depends on other processing variables and the formulation that the manufacturer uses, since wheat flour is the significant and most relevant ingredient used for manufacturing flour tortillas, its quality dramatically affects the final product quality.
- Considering the results of the current available literature, in addition to flour protein quantity and quality, damage starch content should be included to tortilla flour quality specifications.





Increased protein content and dough strength significantly decreased tortilla diameter, but improved tortilla shelf life. Medium protein content and dough strength were ideal to produce good quality tortillas. Four wheat lines (among 131 wheat lines from 1995 to 2007 in SRPN) showed premium tortilla quality, which could be important parental lines in a tortilla wheat breeding program (Zhang et al., 2021).

Eastern US soft winter (ESW) wheat, with high protein content and strength, appears to be suitable for making tortillas. The lower damaged starch content and polyphenol oxidase activity in ESW wheat when compared with hard wheat flour would be advantageous because ESW wheat would require less water for dough preparation and would produce larger and brighter-colored tortillas (Ma and Baik, 2023).







The future of tortilla industry is bright because the plant breeding of high producing nutritionally enhanced both wheat and corn genotypes, novel ingredients, new processing technologies as well as enhanced production equipment and quality testing instruments!







### **Our Mission**

### We craft assurance.

## For our partners. For their customers.

KPM Analytics was formed in 2015 to bring together a portfolio of analytical instrument companies focused on the **food**, **agriculture**, and **environmental** sectors.

Currently, KPM Analytics consists of 9 companies and continues to expand its portfolio of industry experts and premium brands with back-to-back acquisitions.

#### **PURPOSE**

We provide premium quality assurance equipment to food producers through expert craftsmanship and intimate knowledge of their business needs.

#### VISION

Food and agricultural brands the world over will grow stronger from our driven, dedicated, and caring approach to solving their challenges, enabling KPM to become the global industry leader.


#### MISSION

To provide the best solutions for helping our partners control their product quality, scale capacity, and protect their brands.



## Timeline and Progression of KPM Analytics

#### KPM Analytics brands have long and successful histories servicing customers worldwide





Acquires Process Sensors

Founded in 1996

Based in Milford, MA

On-line moisture gauges


NIR technology

### Industries We Serve

KPM's main focus is serving food producers, helping them ensure quality and protect their brand.

KPM is leading the industry for quality solutions at all stages of production.







Our product lines are also widely used in agriculture and feed and forage to measure critical quality parameters.

Environmental, chemistry and industrial industries benefit from our accurate lab and sensing technologies.





## A Global Team

KPM customers are supported by our global sales, service and authorized distribution network.





## Solutions for Every Stage of Production

